

## New Products for LC/MS Applications

- Mobile Phase Blends, OPTIMA LC/MS Grade
- Solvents, OPTIMA LC/MS Grade:  
Acetonitrile, Methanol, 2-Propanol (IPA) and Water
- Reagents, OPTIMA LC/MS Grade:  
Formic Acid, Trifluoroacetic Acid, Ammonium Acetate,  
Ammonium Formate and Acetic Acid



# Mobile Phase Blends, OPTIMA® LC/MS Grade

Fisher Chemical OPTIMA LC/MS solvents have set the standard of excellence for consistent, reproducible performance in the mobile phase of LC/MS. Now, these same high purity solvents are pre-blended with OPTIMA LC/MS modifiers such as formic acid (FA) or trifluoroacetic acid (TFA) to provide ready-to-use aqueous and organic mobile phase blends for LC/MS and LC/UV applications.

In recent laboratory experiments, Fisher Chemical OPTIMA LC/MS solvent blends were shown to meet the stringent purity requirements of LC/MS by providing a consistent concentration of FA or TFA, a very low mass baseline (noise level), exceptionally low metal ion content, and very low LC/UV background. Moreover, the protease-free specification of the aqueous blends is important for proteomics research since peptides/proteins could be degraded if the mobile phase solution is contaminated with protease.

Fisher Chemical OPTIMA LC/MS pre-blended solvents are manufactured in facilities with an ISO 9001:2008 certified quality system to ensure optimum quality and product uniformity.



## KEY FEATURES

- Extensive functional testing to ensure
  - Low mass spectrometry background noise (LC/MS)
  - Minimal metal impurities
  - Lowest impurity background using diode array detection (LC/UV)
- Lot-to-lot consistency
- OPTIMA LC/MS mobile phase extends LC/MS column life due to low impurity levels and low residue value
- Very effective in solubilizing hydrophobic polypeptides
- Only supplier to claim absence of proteases in aqueous blends

## ADVANTAGES

- Ready-to-use
- Avoids the possibility of contamination that can occur with in-house blended solvents
- Eliminates the need to clean glassware or measure corrosive acids
- Eliminates batch-to-batch variation
- Reduces overhead costs associated with preparing blends
- Reduces many of the safety risks associated with storing, blending, and disposing of hazardous solvents and acids
- Avoids human error during preparation of blends

## APPLICATIONS

- Proteomics
- Pharmaceutical Research
- Drug Discovery
- Biomedical Research

OPTIMA LC/MS AQUEOUS BLENDS			
Cat.No.	Product Description	Pack size	Box qty
LS118-500	0.1% FA in Water, OPTIMA LC/MS	500 mL	6 x 500 mL
LS118-1	0.1% FA in Water, OPTIMA LC/MS	1 L	6 x 1 L
LS118-212	0.1% FA in Water, OPTIMA LC/MS	2.5 L	4 x 2.5 L
LS118-4	0.1% FA in Water, OPTIMA LC/MS	4 L	4 x 4 L
LS119-500	0.1% TFA in Water, OPTIMA LC/MS	500 mL	6 x 500 mL
LS119-1	0.1% TFA in Water, OPTIMA LC/MS	1 L	6 x 1 L
LS119-212	0.1% TFA in Water, OPTIMA LC/MS	2.5 L	4 x 2.5 L
LS119-4	0.1% TFA in Water, OPTIMA LC/MS	4 L	4 x 4 L

OPTIMA LC/MS ORGANIC BLENDS			
Cat.No.	Product Description	Pack size	Box qty
LS120-500	0.1% FA in Acetonitrile, OPTIMA LC/MS	500 mL	6 x 500 mL
LS120-1	0.1% FA in Acetonitrile, OPTIMA LC/MS	1 L	6 x 1 L
LS120-212	0.1% FA in Acetonitrile, OPTIMA LC/MS	2.5 L	4 x 2.5 L
LS120-4	0.1% FA in Acetonitrile, OPTIMA LC/MS	4 L	4 x 4 L
LS121-500	0.1% TFA in Acetonitrile, OPTIMA LC/MS	500 mL	6 x 500 mL
LS121-1	0.1% TFA in Acetonitrile, OPTIMA LC/MS	1 L	6 x 1 L
LS121-212	0.1% TFA in Acetonitrile, OPTIMA LC/MS	2.5 L	4 x 2.5 L
LS121-4	0.1% TFA in Acetonitrile, OPTIMA LC/MS	4 L	4 x 4 L

For more information, including a white paper titled 'Optimizing Mobile Phase Solvent Purity for LCMS', go to [www.FisherLCMS.com](http://www.FisherLCMS.com)

# Mobile Phase Blends, OPTIMA LC/MS Grade

## OPTIMA LC/MS Aqueous Blends LS118

0.1% FA IN WATER

OPTIMA LC/MS	
Product Description	OPTIMA LC/MS Water with 0.1% (v/v) Formic Acid

Cat.No. and Pack Size:

LS118-500, 500 mL  
LS118-1, 1 L  
LS118-212, 2.5 L  
LS118-4, 4 L

Assay % (V/V)	0.095 to 0.105%
Color (APHA), Max	10

Optical Abs, wavelength, nm	Au, Max
254	0.01
230	0.55
220	0.85
210	1.25

LC/MS Gradient Suitability	Pass test
Residue After Evaporation (max)	1 ppm
Protease	Not Detected

Ionic Impurities	ppb, Max
Aluminum (Al)	20 ppb
Calcium (Ca)	50 ppb
Copper (Cu)	10 ppb
Iron (Fe)	10 ppb
Lead (Pb)	10 ppb
Magnesium (Mg)	10 ppb
Manganese (Mn)	10 ppb
Nickel (Ni)	10 ppb
Potassium (K)	20 ppb
Silver (Ag)	10 ppb
Sodium (Na)	50 ppb
Zinc (Zn)	20 ppb

## OPTIMA LC/MS Aqueous Blends LS119

0.1% TFA IN WATER

OPTIMA LC/MS	
Product Description	OPTIMA LC/MS Water with 0.1% (v/v) Trifluoroacetic Acid

Cat.No. and Pack Size:

LS119-500, 500 mL  
LS119-1, 1 L  
LS119-212, 2.5 L  
LS119-4, 4 L

Assay % (V/V)	0.095 to 0.105%
Color (APHA), Max	10

Optical Abs, wavelength, nm	Au, Max
254	0.003
230	0.06
220	0.18
210	0.55

LC/MS Gradient Suitability	Pass test
Residue After Evaporation (max)	1 ppm
Protease	Not Detected

Ionic Impurities	ppb, Max
Aluminum (Al)	20 ppb
Calcium (Ca)	50 ppb
Copper (Cu)	10 ppb
Iron (Fe)	10 ppb
Lead (Pb)	10 ppb
Magnesium (Mg)	10 ppb
Manganese (Mn)	10 ppb
Nickel (Ni)	10 ppb
Potassium (K)	20 ppb
Silver (Ag)	10 ppb
Sodium (Na)	50 ppb
Zinc (Zn)	20 ppb

## OPTIMA LC/MS Organic Blends LS120

0.1% FA IN ACETONITRILE

OPTIMA LC/MS	
Product Description	OPTIMA LC/MS Acetonitrile with 0.1% (v/v) Formic Acid

Cat.No. and Pack Size:

LS120-500, 500 mL  
LS120-1, 1 L  
LS120-212, 2.5 L  
LS120-4, 4 L

Assay % (V/V)	0.095 to 0.105%
Color (APHA), Max	10

Optical Abs, wavelength, nm	Au, Max
254	0.03
230	0.75
220	1.25
210	1.3

LC/MS Gradient Suitability	Pass test
Residue After Evaporation (max)	1 ppm
Water (max)	0.01%

Ionic Impurities	ppb, Max
Aluminum (Al)	25 ppb
Calcium (Ca)	50 ppb
Copper (Cu)	10 ppb
Iron (Fe)	10 ppb
Lead (Pb)	10 ppb
Magnesium (Mg)	10 ppb
Manganese (Mn)	10 ppb
Nickel (Ni)	10 ppb
Potassium (K)	20 ppb
Silver (Ag)	10 ppb
Sodium (Na)	50 ppb
Zinc (Zn)	20 ppb

## OPTIMA LC/MS Organic Blends LS121

0.1% TFA IN ACETONITRILE

OPTIMA LC/MS	
Product Description	Optima LC/MS Acetonitrile with 0.1% (v/v) Trifluoroacetic Acid

Cat.No. and Pack Size:

LS121-500, 500 mL  
LS121-1, 1 L  
LS121-212, 2.5 L  
LS121-4, 4 L

Assay % (V/V)	0.095 to 0.105%
Color (APHA), Max	10

Optical Abs, wavelength, nm	Au, Max
254	0.03
230	0.4
220	0.55
210	0.6

LC/MS Gradient Suitability	Pass test
Residue After Evaporation (max)	1 ppm
Water (max)	0.01%

Ionic Impurities	ppb, Max
Aluminum (Al)	25 ppb
Calcium (Ca)	50 ppb
Copper (Cu)	10 ppb
Iron (Fe)	10 ppb
Lead (Pb)	10 ppb
Magnesium (Mg)	10 ppb
Manganese (Mn)	10 ppb
Nickel (Ni)	10 ppb
Potassium (K)	20 ppb
Silver (Ag)	10 ppb
Sodium (Na)	50 ppb
Zinc (Zn)	20 ppb

## Catalogue Cross-References to Competitor Products

Water with 0.1% (v/v) Formic Acid				
<b>Brand</b>	500 mL	1 L	2.5 L	4 L
Fisher Chemical	LS118-500	LS118-1	LS118-212	LS118-4
Fluka	N/A	N/A	34673-2.5L-R	N/A
JT Baker	N/A	9834-02	N/A	9834-03
Biosolve	N/A	232441, 1 L	232441, 2.5 L	N/A
Water with 0.1% (v/v) Trifluoroacetic Acid				
<b>Brand</b>	500 mL	1 L	2.5 L	4 L
Fisher Chemical	LS119-500	LS119-1	LS119-212	LS119-4
Fluka	N/A	N/A	34978-2.5L-R	N/A
JT Baker	N/A	9836-02	N/A	9836-03
Biosolve	N/A	232741	N/A	N/A

Acetonitrile with 0.1% (v/v) Formic Acid				
<b>Brand</b>	500 mL	1 L	2.5 L	4 L
Fisher Chemical	LS120-500	LS120-1	LS120-212	LS120-4
Fluka	N/A	N/A	34668-2.5L-R	N/A
JT Baker	N/A	9832-02	N/A	9832-03
Biosolve	N/A	019341, 1 L	019341, 2.5 L	N/A
Acetonitrile with 0.1% (v/v) Trifluoroacetic Acid				
<b>Brand</b>	500 mL	1 L	2.5 L	4 L
Fisher Chemical	LS121-500	LS121-1	LS121-212	LS121-4
Fluka	N/A	N/A	34976-2.5L-R	N/A
JT Baker	N/A	9835-02	N/A	9835-03
Biosolve	N/A	019541, 1 L	019541, 2.5 L	N/A

# Solvents, OPTIMA LC/MS Grade

The joining of liquid chromatography (LC) with mass spectrometry (MS) has become an indispensable tool for various fields of research. The value of LC/MS derives from its ability to combine separation chemistry with selective mass ion detection. As instrumentation advances lead to ever-lower analyte detection limits, it is crucial for the chromatographer to consider the level of purity when selecting appropriate solvents for use in the LC/MS mobile phase. The Fisher Chemical product line offers superior high-purity solvents designed to meet the required purity level of advanced LC/MS systems.



## OPTIMA LC/MS Acetonitrile A955

Assay (by GC), min	99.9%
--------------------	-------

Optical Abs, wavelength, nm	Au, Max
280	0.005
254	0.005
230	0.01
225	0.015
220	0.015
215	0.025
210	0.03
205	0.04
200	0.05
195	0.15
190	1.00

LC/UV Gradient Suitability, nm	Single max peak (Au)
254	0.0005
210	0.002

LC/MS at Positive Mode Max	50 ppb Reserpine
LC/MS at Negative Mode Max	50 ppb Aldicarb
Water (KF) %	0.01
Residue after Evaporation, ppm Max	0.8

Trace Ionic Impurities	ppb, Max
Aluminum (Al)	25
Barium (Ba)	5
Cadmium (Cd)	5
Calcium (Ca)	25
Chromium (Cr)	5
Cobalt (Co)	5
Copper (Cu)	5
Iron (Fe)	5
Lead (Pb)	5
Manganese (Mn)	5
Magnesium (Mg)	10
Nickel (Ni)	5
Potassium (K)	10
Silver (Ag)	5
Sodium (Na)	50
Tin (Sn)	5
Zinc (Zn)	10

Titrate Acid, mEq/g	0.008
Titrate Base, mEq/g	0.0006

## OPTIMA LC/MS Methanol A456

Assay (by GC), min	99.9%
--------------------	-------

Optical Abs, wavelength, nm	Au, Max
280	0.005
260	0.005
254	0.01
230	0.1
220	0.2
214	0.4
210	0.5

LC/UV Gradient Suitability, nm	Single max peak (Au)
254	0.001
220	0.005

LC/MS at Positive Mode Max	50 ppb Reserpine
LC/MS at Negative Mode Max	50 ppb Aldicarb
Water (KF) %	0.02
Residue after Evaporation, ppm, Max	1

Trace Ionic Impurities	ppb, Max
Aluminum (Al)	10
Barium (Ba)	10
Cadmium (Cd)	10
Calcium (Ca)	20
Chromium (Cr)	10
Cobalt (Co)	10
Copper (Cu)	10
Iron (Fe)	10
Lead (Pb)	10
Magnesium (Mg)	10
Manganese (Mn)	10
Nickel (Ni)	10
Potassium (K)	10
Silver (Ag)	10
Sodium (Na)	50
Tin (Sn)	10
Zinc (Zn)	10

Titrate Acid, mEq/g	0.0003
Titrate Base, mEq/g	0.0002

## OPTIMA LC/MS Water W6

Optical Abs, wavelength, nm	Au, Max
280	0.005
260	0.005
254	0.005
240	0.01
230	0.01
220	0.01
210	0.01

LC/UV Gradient Suitability, nm	Single max peak (Au)
254	0.0005
210	0.005

Protease	Not Detected
LC/MS at Positive Mode Max	50 ppb Reserpine
LC/MS at Negative Mode Max	50 ppb Aldicarb
Residue after Evaporation, ppm, Max	1

Trace Ionic Impurities	ppb, Max
Aluminum (Al)	10
Barium (Ba)	10
Cadmium (Cd)	10
Calcium (Ca)	20
Chromium (Cr)	10
Cobalt (Co)	10
Copper (Cu)	10
Iron (Fe)	10
Lead (Pb)	10
Magnesium (Mg)	10
Manganese (Mn)	10
Nickel (Ni)	10
Potassium (K)	10
Silver (Ag)	10
Sodium (Na)	20
Tin (Sn)	10
Zinc (Zn)	10
Total Halogens (as chloride)	Not Detected

## OPTIMA LC/MS 2-Propanol A461

Assay	99.9%
-------	-------

Color (APHA), Max	5
-------------------	---

Optical Abs, wavelength, nm	Au, Max
254	0.005
230	0.05
220	0.1
210	0.4

LC/MS Suitability	Pass Test
-------------------	-----------

Water (KF) %	0.05
--------------	------

Residue after Evaporation, ppm, Max	1
-------------------------------------	---

Titrate Acid or Base meq/g	0.0001
----------------------------	--------

Trace Ionic Impurities	ppb, Max
Aluminum (Al)	10
Calcium (Ca)	10
Copper (Cu)	5
Iron (Fe)	5
Lead (Pb)	5
Magnesium (Mg)	5
Manganese (Mn)	5
Nickel (Ni)	5
Potassium (K)	10
Silver (Ag)	5
Sodium (Na)	50
Zinc (Zn)	10

## Catalogue Cross-References to Competitor Products

Acetonitrile					
Brand	500 mL	1 L	2.5 L	4 L	
Fisher Chemical	A955-500	A955-1	A955-212	A955-4	
EMD	NA	AX0156-6	NA	AX0156-1	
Merck	NA	1000291000	1000292500	NA	
JT Baker EU	NA	98531000	NA	NA	
JT Baker US	NA	9829-2	NA	9829-3	
Fluka	NA	34967-1L	34967-2.5L	34967-4X4L	
BioSolve	12041	12041	NA	NA	
Methanol					
Brand	500 mL	1 L	2.5 L	4 L	
Fisher Chemical	A456-500	A456-1	A456-212	A456-4	
EMD	NA	MX0486-6	NA	MX0496-6	
Merck	NA	1060351000	1060352500	NA	
JT Baker	NA	98631000	NA	NA	
		9830-2		9830-3	
Fluka	NA	34966-1L	34966-2.5L	34966-4X4L	
BioSolve	136841	136841	NA	NA	

2-Propanol					
Brand	500 mL	1 L	2.5 L	4 L	
Fisher Chemical	A461-500	A461-1	A461-212	A461-4	
EMD	NA	NA	NA	NA	
Merck	NA	NA	NA	NA	
JT Baker EU	NA	98301000	NA	NA	
JT Baker US	NA	9627-2	NA	9627-3	
Fluka	NA	34965-1L	34965-2.5L	NA	
BioSolve	162641	162641	162641	NA	
Water					
Brand	500 mL	1 L	2.5 L	4 L	
Fisher Chemical	W6-500	W6-1	W6-212	W6-4	
EMD	NA	WX0001-6	NA	WX0001-1	
Merck	NA	NA	NA	NA	
JT Baker	NA	98231000	NA	NA	
		9831-2		9831-3	
Fluka	NA	39253-1L	NA	39253-4X4L	
BioSolve	NA	232141	232141	NA	

# Reagents, OPTIMA LC/MS Grade

The Fisher Chemical OPTIMA LC/MS product line now includes neat additives such as A116 Trifluoroacetic Acid (TFA) and A117 Formic Acid (FA), coupled with their ready-to-use aqueous and organic blends. These mobile phase additives are use-tested, and our data will show that their mass baselines are very low for positive mode in TIC with single quadrupole and ion trap mass detector, very low LC-UV response using diode array detector, and exceptionally low metal ion content with ICP-MS which makes MS interpretation easier due to reduced metal adduct formation. Additionally, for proteomics research the neat additives and blends are protease free which ensures that target peptides or proteins are not destroyed by these pervasive hydrolyzing enzymes.

When compared to the competition, Fisher Chemical OPTIMA LC/MS Trifluoroacetic Acid (A116), Formic Acid (A117), and their associated mobile phase blends exhibit:

- Low mass spectrometry background noise
- Lowest impurity background using diode array detection (LC-UV)
- Minimal metal impurities
- Aqueous mobile phase blends free of proteases
- Lot-to-lot consistency
- Extended column life due to low impurity levels using OPTIMA LC/MS solvents and additives



## Formic Acid, OPTIMA LC/MS Grade

Fisher Chemical OPTIMA LC/MS Grade Formic Acid is an ultrapure reagent used as an additive for the formulation of solvent blends for the mobile phase in LC/MS applications. Formic acid is often added to the LC/MS mobile phase due to its ability to reduce MS signal suppression and thus enhance the signal detection limit.

For preparing mobile phase blends, A117-50 Formic Acid can be used with OPTIMA LC/MS Solvents such as A955 Acetonitrile, W6 Water, and A456 Methanol.



PRODUCT SPECIFICATIONS	A117
	FORMIC ACID
<b>Appearance</b>	<b>Clear, colorless, fuming liquid</b>
Assay	99.5% min
Color(APHA), Max	10

Optical Absorbance (0.1% aqueous), au, max	
254 nm	0.01
230 nm	0.55
220 nm	0.85
210 nm	1.25

LCMS Suitability (at 0.1% aqueous)	Pass Test
------------------------------------	-----------

Trace Ionic Impurities (at 0.1% aqueous), ppb, max	
Aluminum (Al)	20
Calcium (Ca)	50
Copper (Cu)	10
Iron (Fe)	10
Lead (Pb)	10
Magnesium (Mg)	10
Manganese (Mn)	10
Nickel(Ni)	10
Potassium (K)	20
Silver (Ag)	10
Sodium (Na)	50
Zinc (Zn)	20

Residue after evaporation (at 0.1% aqueous)	1 ppm max
Chloride, max	5 ppm
Sulfate (SO <sub>4</sub> ), max	50 ppm
Sulfite (SO <sub>3</sub> ), max	20 ppm
Water (KF), max	

## PACKAGING INNOVATIONS

- Safety first - packaged in HDPE bottle to avoid possible breakage of glass bottle due to pressure buildup from carbon monoxide which is a natural decomposition product of formic acid
- Proprietary surface treatment applied to HDPE bottle to create a barrier between bottle and formic acid thus preventing contamination by plasticizers

Cat.No.	Product Description	Pack Size
<b>A117-50</b>	Formic Acid, Optima LC/MS	50 mL

For smaller volumes we now offer Formic Acid in ampules		
Cat.No.	Product Description	Pack Size
<b>A117-10X1AMP</b>	Formic Acid, OPTIMA LC/MS	10 x 1 mL
<b>A117-1AMP</b>	Formic Acid, OPTIMA LC/MS	1 mL
<b>A117-05AMP</b>	Formic Acid, OPTIMA LC/MS	0.5 mL
<b>A117-2AMP</b>	Formic Acid, OPTIMA LC/MS	2 mL



# Reagents, OPTIMA LC/MS Grade

## Trifluoroacetic Acid, OPTIMA LC/MS Grade

Fisher Chemical OPTIMA LC/MS Grade Trifluoroacetic Acid is an ultrapure reagent used as an additive for the formulation of solvent blends for the mobile phase in LC/MS applications.

Trifluoroacetic Acid is often added to the mobile phase to enhance chromatographic separations. When reversed-phase chromatography is applied to protein and peptide analysis, TFA is frequently used as an ion-pairing reagent in mobile phases. Although TFA is known to suppress MS ionization, use of a low concentration of TFA can enhance retention and improve peak shape of a product during the analytical separations. For preparing mobile phase blends, A116-50 Trifluoroacetic Acid can be used with OPTIMA LC/MS Solvents such as A955 Acetonitrile, W6 Water, and A456 Methanol.



### TFA AND FA KEY FEATURES

- Also available in convenient 0.5 mL, 1 mL and 2 mL ampules for easy dilutions
- Ampules are prescored for easy opening - no need to file the ampule
- Minimal metal impurities
- Low mass spectrometry background noise
- Lowest impurity background using diode array detection (LC-UV)
- Extended column life due to low impurity levels using OPTIMA LC/MS mobile phase prepared with A116 Trifluoroacetic Acid and A117 Formic Acid

Cat.No.	Product Description	Pack Size
<b>A116-50</b>	Trifluoroacetic Acid, Optima LC/MS	50 mL

For smaller volumes we now offer Trifluoroacetic Acid in ampules		
Cat.No.	Product Description	Pack Size
<b>A116-10X1AMP</b>	Trifluoroacetic Acid, OPTIMA LC/MS	10 x 1 mL
<b>A116-1AMP</b>	Trifluoroacetic Acid, OPTIMA LC/MS	1 mL
<b>A116-05AMP</b>	Trifluoroacetic Acid, OPTIMA LC/MS	0.5 mL
<b>A116-2AMP</b>	Trifluoroacetic Acid, OPTIMA LC/MS	2 mL

PRODUCT SPECIFICATIONS	A116
	TRIFLUOROACETIC ACID
<b>Appearance</b>	<b>Clear, colorless, fuming liquid</b>
Assay	99.5% min
Color(APHA), Max	10

Optical Absorbance (as is), au, max	
300 nm	0.03
275 nm	0.06

Optical Absorbance (0.1% aqueous), au, max	
254 nm	0.003
230 nm	0.06
220 nm	0.18
210 nm	0.55

LCMS Suitability (at 0.1% aqueous)	Pass Test
------------------------------------	-----------

Trace Ionic Impurities (at 0.1% aqueous), ppb, max	
Aluminum (Al)	20
Calcium (Ca)	50
Copper (Cu)	10
Iron (Fe)	10
Lead (Pb)	10
Magnesium (Mg)	10
Manganese (Mn)	10
Nickel(Ni)	10
Potassium (K)	20
Silver (Ag)	10
Sodium (Na)	50
Zinc (Zn)	20

Residue after evaporation (at 0.1% aqueous)	1 ppm max
Chloride, max	
Sulfate (SO4), max	
Sulfite (SO3), max	
Water (KF), max	0.05%

# Reagents, OPTIMA LC/MS Grade

## New! Acetic Acid, Ammonium Acetate, and Ammonium Formate

Achieving the full analytical potential of LC/MS instruments requires the use of highly purified solvents and ionic additives. Fisher Chemical OPTIMA LC/MS Ammonium Formate, Ammonium Acetate, and Acetic Acid are ultrapure reagents used to modify the LC/MS mobile phase in order to improve chromatographic peak shape and provide stable analyte signals in the MS detector.



Ammonium acetate and ammonium formate are volatile salts often used for improving ionization under neutral conditions in either the positive or negative ESI-MS mode. Acetic acid is commonly used for enhancing chromatography separation under acidic conditions.

For preparing mobile phase blends, A113-50 Acetic Acid, A114-50 Ammonium Acetate, and A115-50 Ammonium Formate can be added to Fisher Chemical OPTIMA LC/MS solvents such as A955 Acetonitrile, W6 Water, and A456 Methanol.

CATALOG NO.	PRODUCT DESCRIPTION	PACK SIZE	PACKAGING
A113-50	Acetic Acid, OPTIMA LC/MS	50 mL	Treated HDPE
A114-50	Ammonium Acetate, OPTIMA LC/MS	50 g	Amber Glass
A115-50	Ammonium Formate, OPTIMA LC/MS	50 g	Amber Glass

PRODUCT SPECIFICATIONS	A113	A114	A115
	ACETIC ACID	AMMONIUM ACETATE	AMMONIUM FORMATE
Appearance	Clear, colorless fuming liquid	White fine crystals w/ lumps	White fine crystals w/ lumps
Assay, % min	99.7% min	99% min	99% min
Color (APHA), max	10 max		

Optical Absorbance, au, max			
350 nm	0.02 au max		
280 nm	0.05 au max		
255 nm	1.0 au max		

LC/MS Gradient Suitability (at 0.1% aqueous)	A113	A114	A115
	Pass test	Pass test	Pass test

Trace Ionic Impurities (at 0.1% aqueous), ppb, max			
Aluminum (Al)	20 ppb	20 ppb	20 ppb
Calcium (Ca)	50 ppb	50 ppb	50 ppb
Copper (Cu)	10 ppb	10 ppb	10 ppb
Iron (Fe)	10 ppb	10 ppb	10 ppb
Lead (Pb)	10 ppb	10 ppb	10 ppb
Magnesium (Mg)	10 ppb	10 ppb	10 ppb
Manganese (Mn)	10 ppb	10 ppb	10 ppb
Nickel (Ni)	10 ppb	10 ppb	10 ppb
Potassium (K)	20 ppb	20 ppb	20 ppb
Silver (Ag)	10 ppb	10 ppb	10 ppb
Sodium (Na)	50 ppb	50 ppb	50 ppb
Zinc (Zn)	20 ppb	20 ppb	20 ppb

Residue after Evaporation, ppm, max	5 ppm max		
Residue on Ignition, % max		0.01% max	0.005% max
Chloride, % max	0.00004% max	0.0005% max	0.001% max
Nitrate (NO <sub>3</sub> ), % max		0.001% max	
Sulfate (SO <sub>4</sub> ), % max	0.00004% max		0.005% max
Protease (at 0.1% aqueous)	Not detected	Not detected	Not detected
Water (KF), % max	0.1% max	1% max	2% max

Shelf Life	A113	A114	A115
	1 year	2 years	3 years

## KEY FEATURES

- Minimal metal impurities
- Low mass spectrometry background noise
- Lowest impurity background using diode array detection (LC/UV)
- Lot-to-lot consistency
- Extended column life due to low impurity levels using OPTIMA LC/MS mobile phase solvents and additives

## PACKAGING INNOVATION

- Proprietary surface treatment applied to HDPE bottle to create a barrier between bottle and acetic acid, thus preventing contamination by plasticizers

## New! OPTIMA LC/MS Acetic Acid is available in convenient glass ampules

Fisher Chemical OPTIMA LC/MS Acetic Acid is now custom-packaged in amber borosilicate ampules under inert atmospheric conditions to provide the freshest additive for preparing aqueous and organic mobile phase blends.



CATALOG NO.	PRODUCT DESCRIPTION	PACK SIZE
A113-1AMP	Acetic Acid, OPTIMA LC/MS	1 mL ampule
A113-10X1AMP	Acetic Acid, OPTIMA LC/MS	10 x 1 mL ampule

## ADVANTAGES

- Convenient - ready-to-use ampules in 1 ml size
- Safe - eliminates the need to measure corrosive acids
- Reproducible - reduces batch-to-batch variation since ampules are filled to +/- 0.05%
- Expedient - 1 ml ampule is used to prepare in a matter of seconds 1 L of fresh 0.1% v/v acetic acid mobile phase blend

**GLOBAL LOCATIONS**
**AMERICAS**
**Canada**

Fisher Scientific Canada  
112 Colonnade Road  
Ottawa, Ontario  
Post Code: K2E 7L6  
Toll-Free Number: 800-234-7437  
Fax: 800-463-2996  
www.fishersci.ca

**Latin America**

Fisher Scientific Global Export, Latin America  
3970 Johns Creek Court  
Suite 500  
Suwanee, GA  
Post Code: 30024  
Toll-Free Number: 770-871-4725  
Fax: 770-871-4726  
www.fishersci.com

**United States**

Fisher HealthCare  
9999 Veterans Memorial Drive  
Houston, TX  
Post Code: 77038  
Toll-Free Number: 800-640-0640  
Fax: 800-290-0290  
www.fisherhealthcare.com

Fisher Scientific  
300 Industry Drive  
Pittsburgh, PA  
Post Code: 15275  
Toll-Free Number: 800-766-7000  
Fax: 800-926-1166  
www.fishersci.com

**ASIA**
**China**

Fisher Scientific Worldwide (Shanghai) Co., Ltd  
Hotline: +86 400 881 5117  
website: www.fishersci.com.cn  
E-mail: sales.china@thermofisher.com

Shanghai office  
2F, Building 7  
No. 27 Xin Jingqiao Rd  
Shanghai 201206  
P.R. China  
Tel: +86 21-52581100; +86 21-61925688  
Fax: +86 21-52580119

Beijing office  
Units 702-715, 7th Floor,  
Tower West, Yonghe Plaza  
No.28 Andingmen East Street  
Beijing 100007  
P.R. China  
Tel: +86 10-84193588  
Fax: +86 10-84193580

Guangzhou office  
Room 2405-2406, Time Estate Center  
No.410-412 Middle DongFeng Road  
Guangzhou 510030  
P.R. China  
Tel: +86 20-83145288  
Fax: +86 20-38771941

**India**

Fisher Scientific (Part of Thermo Fisher Scientific)  
403-404, B-Wing, Delphi  
Hiranandani Business Park  
Powai, Mumbai-400 076  
India  
Toll-Free Number: 1 800 209 7001  
Fax: 91 22 6680 3001/3002  
qfc.customer@thermofisher.com  
www.fishersci.in

**Japan**

Fisher Scientific Japan  
Thermo Fisher Scientific K.K.  
C-2F, 3-9 Moriya-cho  
Kanawaga-ku, Yokohama  
221-0022 Japan  
Tel: 81 45 450 6310  
Fax: 81 45 450 6316  
support@fishersci.co.jp  
www.fishersci.co.jp

**Korea**

Fisher Scientific Korea  
Sambu Bldg. 13F  
676 Yeoksam-dong,  
Kangnam-Gu  
Seoul 135-979, Korea  
Dir (02) 527-0300  
Customer Service Center: (02) 527-0300  
Fax: (02) 527-0311  
sales@fishersci.co.kr  
www.acros.co.kr

**Malaysia**

Fisher Scientific Malaysia Sdn Bhd  
No. 3 Jalan Sepadu 25/123  
Taman Perindustrian Axis Seksyen 25  
40400 Shah Alam  
Selangor Darul Ehsan Malaysia  
Technical Service Hotline: 1-300-88-7868  
Tel: (603) 51218888  
Fax: (603) 51218899  
enquiry.my@thermofisher.com  
www.fishersci.com.my

**Singapore**

Fisher Scientific Singapore  
Fisher Scientific Pte Ltd.  
8 Pandan Crescent  
LL4, #05-05 UE Tech Park  
Singapore 128464  
Tel: (65) 6873 6006  
Fax: (65) 6873 5005  
enquiry.sg@thermofisher.com  
www.fishersci.com.sg

**EUROPE**
**Austria**

Fisher Scientific (Austria) GmbH  
Rudolf von Alt-Platz 1  
A-1030 Wien  
Toll-Free Number: 0800 20 88 40  
Fax: 0800 20 66 90  
info.austria@thermofisher.com  
www.at.fishersci.com

**Belgium**

Fisher Scientific  
BP 567  
B-7500 Tournai 1  
Tel: 056 260 260  
Fax: 056 260 270  
be.fisher@thermofisher.com  
www.be.fishersci.com

**Czech Republic**

Fisher Scientific, spol. s r.o.  
Kosmonautů 324  
Pardubice  
CZ-530 09  
Tel: 466 798 230  
Fax: 466 435 008  
info.cz@thermofisher.com  
www.thermofisher.cz

**Denmark**

Fisher Scientific  
Industrivej 3  
DK-3550 Slangerup  
Tel: +45 70 27 99 20  
Fax: +45 70 27 99 29  
kundeservice@thermofisher.com  
www.fishersci.dk

**Finland**

Fisher Scientific Oy  
Ratastie 2  
FI-01620 Vantaa  
Tel: +358 802 76 280  
Fax: +358 802 76 235  
fisher.fi@thermofisher.com  
www.fishersci.fi

**France**

Fisher Scientific  
Parc d'Innovation BP 50111  
67403 Illkirch Cedex  
Tel: 03 88 67 53 20  
Fax: 03 88 67 11 68  
fr.commande@thermofisher.com  
www.fr.fishersci.com

**Germany**

Fisher Scientific GmbH  
Im Heiligen Feld 17  
D-58239 Schwerte  
Tel: +49 (0)2304 932-5  
Fax: +49 (0)2304 932-950  
info.germany@thermofisher.com  
www.de.fishersci.com

**Ireland**

Fisher Scientific Ireland  
Suite 3 Plaza 212  
Blanchardstown Corporate  
Park 2  
Ballycoolin  
Dublin 15  
Tel: +353 01 885 5854  
Fax: +353 01 899 1855  
fsie.sales@thermofisher.com  
www.ie.fishersci.com

**Italy**

Fisher Scientific  
Tel: 02 953 28 258  
Fax: 02 953 27 374  
it.fisher@thermofisher.com  
www.it.fishersci.com

**The Netherlands**

Fisher Scientific  
Scheepsbouwersweg 1b  
Postbus 4  
1120 AA Landsmeer  
Tel: 020 487 70 00  
Fax: 020 487 70 70  
nl.info@thermofisher.com  
www.nl.fishersci.com

**Norway**

Fisher Scientific  
Frysjavaen 33E  
0884 Oslo  
Tel: +47 22 95 59 59  
Fax: +47 22 95 59 40  
fisher.no@thermofisher.com  
www.fishersci.no

**Portugal**

Fisher Scientific  
Rua Pedro Álvares Cabral, n°24, 3°D  
Edifício Euro - Infatado  
2670-391 Loures  
Tel: +351 21 425 33 50/4  
Fax: +351 21 425 33 51  
pt.fisher@thermofisher.com  
www.pt.fishersci.com

**Spain**

Fisher Scientific  
C/ Luis I, 9  
28031 MADRID  
Tel: 91 380 67 10  
Fax: 91 380 85 02  
es.fisher@thermofisher.com  
www.es.fishersci.com

**Sweden**

Fisher Scientific  
Box 9193  
400 94 Göteborg  
Tel: +46 31 352 32 00  
Fax: +46 31 352 32 50  
fisher.se@thermofisher.com  
www.fishersci.se

**Switzerland**

Fisher Scientific  
Wilstrasse 57  
5610 Wohlen  
Tel: 056 618 41 11  
Fax: 056 618 41 41  
info.ch@thermofisher.com  
www.ch.fishersci.com

**United Kingdom**

Fisher Scientific UK  
Bishop Meadow Road  
Loughborough  
Leicestershire LE11 5RG  
Tel: +44 (0)1509 555500  
Fax: +44 (0)1509 555111  
fsuk.sales@thermofisher.com  
www.fisher.co.uk

**Rest of Europe**

Thermo Fisher Scientific  
ENA 23, Zone 1, nr 1350  
Janssen Pharmaceuticaaan 3a  
2440 Geel – Belgium  
Tel: +32 14 57 52 11  
Fax: +32 14 59 26 10  
www.acros.com

**MIDDLE EAST AND AFRICA**
**Fisher Scientific Global Export**

**MEA Office**  
DAFZA  
Office Block 4A #504  
P.O. Box 54710  
Dubai, UAE  
Tel: +971 4609 1594  
Fax: +971 4204 5564  
www.fishersci.com

**OCEANIA**
**Australia**

Thermo Fisher Scientific Australia Pty Ltd  
5 Caribbean Drive  
Scoresby, VIC 3179  
Tel: +1300 735 292  
Fax: +61 3 9763 1169  
AUinfo@thermofisher.com  
www.thermofisher.com.au

**New Zealand**

Thermo Fisher Scientific New Zealand Ltd  
244 Bush Road,  
Albany, Auckland 0632  
Toll-Free Number: 0800 933 966  
Fax: +64 9 980 6788  
infonz@thermofisher.com  
www.thermofisher.co.nz

©2011 Thermo Fisher Scientific Inc.  
All rights reserved. All brands and trademarks  
are part of Thermo Fisher Scientific Inc. and  
its subsidiaries.

BN0501093